



Datasheet

Model No.: MNXT2S-Xxx211214





#### General features





Illustration of single particle sensitivity of Timepix2 device. The tracks of different particles of radiation background (mostly muons and few protons) were recorded in 5 minutes on board of airplane. No noise (clean zero) is seen in dark regions.

The MINIPIX TPX2 is a miniaturized and low-power radiation camera solution that incorporates a single Timepix2 detector with a sensor of customer preference (typically 300 µm thick silicon). The detector features 256 x 256 pixels with a pitch of 55 µm and is capable of single particle counting or particle tracking. The MINIPIX TPX2 utilizes a USB 2.0 interface, allowing for reading of up to 99 frames per second. The energy-sensitive Timepix2 detector brings a new dimension to radiographic images and now also features a new measurement modality - adaptive gain. Adaptive gain helps to improve performance in high-intensity use cases, increasing the dynamic range of the device, and making it an even more versatile and powerful tool for radiation detection.

The MINIPIX TPX2 device is controlled via a USB interface and is compatible with major operating systems such as MS Windows, Mac OS, and LINUX. The system includes free software, PIXET PRO, for detector operation, offering comprehensive functionality and ease of use. With its miniaturized size, low power consumption, and advanced Timepix2 detector technology, the MINIPIX TPX2 is an efficient and effective solution for various radiation detection applications (imaging, XRD, particle tracking etc.)

# Main Features

•	Readout chip type	Timeniy?
		·
•	Pixel size	. 55 x 55 μm <sup>*</sup>
•	Sensor resolution	. 256 x 256 pixels
•	Counter bit depth	. 10/ 14/ 18 bit
•	Dark current	. none
•	Interface	. USB 2.0 (Full-Speed)
•	Maximum frame rate	. 99 fps
•	Dimensions	. 88.9 x 21 x 10 mm
•	Weight	. 30 g

- $^*$  55 x 110  $\mu m$  at the edges and 110 x 110  $\mu m$  at the corners
- <sup>1</sup> MINIPIX TPX2 is not certified dosimetric device. It serves as the first level indicator and monitor of radiation fields allowing identification of a radiation type. Radiation protection of people cannot be based on measurements of MINIPIX TPX2.
- <sup>2</sup> Dynamic range of final picture is theoretically unlimited; the only limiting factor is exposure time.





# Device parameters

## Operating conditions

Symbol	Parameter	Value	Units	Comment
TA	Ambient Temperature Range	0-50	°C	
Ф	Humidity	<80	%	Not condensing
	Altitude	<2000	m	Above sea level
IP	IP rating	IP40		With cover

#### **Electrical Specification**

 $T_A = 25$ °C, USB voltage  $V_{CC} = 4.8V$ 

Symbol	Parameter	Min	Тур	Max	Units	Comment
V <sub>CC</sub>	Supply Voltage	4.0	5.0	5.5	V	
Icc2	Chip active		500	700	mA	
P1	Power Dissipation		2.5	3.5	W	
Bias Voltage Source for Sensor Diode						
V <sub>BIAS</sub>	Bias Voltage	3		200	V	

#### Performance characteristics

Symbol	Parameter	Min	Тур	Max	Units	Comment
f	Frame-rate			99	fps	with USB 2.0 Host
T <sub>READ</sub>	Frame Readout Time <sup>3</sup>		19		ms	




<sup>&</sup>lt;sup>3</sup> During Readout time (or Dead time), no charge is collected from the sensor.



# Sensor parameters

 $T_A = 25^{\circ}C$ 

Symbol	Parameter	Si		Units	Comment	
	Thickness	100	300	500	μm	
σ	Energy resolution of energy discrimination threshold (σ @ 8 keV) 0.36		keV			
σ	Energy resolution in full spectral mode (σ @ 8 keV)	0.6 keV				
σ	Energy resolution in full spectral mode (σ @ 23 keV)	0.9 keV		keV		
σ	Energy resolution in full spectral mode (σ @ 60 keV) 1.4		keV			
	Typical detectable energy range for X-rays 5 to 60		keV	See chart below		
	Pixel size		55		μm	







## Modes and types of readout chip operation

The detector is frame-based: The data from all the pixels are read out after the acquisition time is over.

#### Modalities:

Integral measurement:

During the acquisition, recorded data is integrated and outputted as a single frame.

First hit measurement:

This mode disregards events that take place in the same pixel during the acquisition time, in order to minimize pile-ups.

Counter bit Depth:

Different counter-bit depths can be chosen for certain measurement modes. This enables tailoring the performance for higher frame rates, or better resolution.

Combinations of operation modes and measurement modalities (default cases are highlighted):

Mode	Counter Depth	Energy measurement	Frame Rate	
Counts	14 bits	N/A	64 fps	
Counts	10 bits - high frame rate	11//	99 fps	
Energy		Integrated energy	- 65 fps	
Lifergy	14 bits	1st hit measurement		
Time	14 bits	N/A	65 fps	
Time	10 bits - high frame rate	IV/A	98 fps	
Counts +	10 hits (Enorgy) / 4 hits	Integrated energy	- 61 fps	
Energy	10 bits (Energy) / 4 bits (Counts)	1st hit measurement		
	14 bits (Energy) / 14 bits (Time)	Integrated energy	32 fps	
Energy +	14 orts (Energy)/ 14 orts (Time)	1st hit measurement	32 ip3	
Time	10 bits (Energy) / 18 bits (Time)	Integrated energy	34 fps	
	10 5.03 (Energy) / 10 5.03 (Time)	1 <sup>st</sup> hit measurement	34 ips	






# Vacuum Operation

Advacam detectors are vacuum compatible out of the box. Operate only with air pressure lower than  $10^{-3}$ Pa. Intended for dust-free indoor use.

Make sure to disconnect the device from power during pumping down or venting the vacuum chamber!

## **Device description**



#### **USB** connector

USB type Micro-B, Standard USB 2.0 High-Speed.


The USB cable length should be less than 2m!

For longer connections, a repeater or active cable is suggested.





#### Mechanical dimensions



All dimensions are in mm.

Extreme care must be taken when removing protecting cover and handling the **MINIPIX**<sub>TPX2</sub> without the protecting cover. The warranty does not apply to mechanical damage of the sensor and wirebonds.



<sup>\*</sup> Sensitive surface distance from top of the box is for 300  $\mu m$  sensor thickness.



# **Model Number Codes**

Example:	MNX	T2S	- <u>X</u>	P	3	210520
Device name:						
MNX – MiniPIX						
Device modification:						
T2S – Timepix Standard						
Sensor type:						
P – Planar silicon						
Sensor thickness:						
1 – 100 μm						
3 – 300 μm						
5 – 500 μm						
Device version date:						
YY MM DD						

# Release history

Date	Changes
25/04/2023	New version





# Warning

### Do not touch sensor surface!

#### Instructions for safe use

To avoid malfunction or damage to your MINIPIXTPX2 please observe the following:

- Do not expose to water or moisture.
- Do not disassemble. Wire-bonding connection may be irreversibly damaged.
- Do not insert any object into the sensor window.
- Maximum USB cable length is 2m
- The protection provided by this product may be impaired if it is used in a manner not described in this document

#### Disposal:



Do not dispose these instruments as unsorted municipal waste. Please use separate collection facility to contact the supplier from which the instrument was purchased. Please make sure discarded electrical waste is properly recycled to reduce environment impact

## Copyright

ADVACAM s.r.o. Tel: +420 603 444 112

U Pergamenky 1145/12 Email:

170 00 Praha info@advacam.com

Czech Republic www.advacam.com

